



## ACID STORAGE CABINETS

NuAire Acid Storage Cabinets are constructed from all stress-relieved, fully seam-welded, white corrosive-resistant polypropylene.

The outer cabinet shelf, kick plate and adjustable shelves with drip-proof lips are constructed from ½ inch polypropylene.

Each cabinet incorporates easy-to-clean shelves, handles, hinges, and screws also made of 100% polypropylene.

Adjustable slide front air vents and removable top air vent plugs assure acid storage cabinets have adequate ventilation for safety.

DOOR LOCKS AND CASTERS ARE OPTIONAL.

## ACID STORAGE CABINETS



| MODEL NUMBER | SIZE (W x H x D) |
|--------------|------------------|
| NU-53-3660   | 36" x 60" x 18"  |
| NU-53-4260   | 42" x 60" x 18"  |
| NU-53-4860   | 48" x 60" x 18"  |
| NU-53-3684   | 36" x 84" x 22"  |
| NU-53-4284   | 42" x 84" x 22"  |
| NU-53-4884   | 48" x 84" x 22"  |

**OPTIONS:** Casters, Vent Location, Key Locks

## POLYPROPYLENE TABLES



| MODEL NUMBER | SIZE (W x H x D) |
|--------------|------------------|
| NU-98-3636   | 36" x 36" x 24"  |
| NU-98-4836   | 48" x 36" x 24"  |
| NU-98-6036   | 60" x 36" x 24"  |
| NU-98-7236   | 72" x 36" x 24"  |

**OPTIONS:** Casters, Epoxy Resin Tabletop

## PEG BOARDS



**WHITE/BLACK W/PEGS & 3" DRAIN TROUGH**


| MODEL NUMBER | SIZE (W x H) |
|--------------|--------------|
| NU-56-3030   | 30" x 30"    |
| NU-56-3036   | 30" x 36"    |

WITH NUAIRE YOU KNOW



[www.nuaire.com](http://www.nuaire.com)

2100 Fernbrook Lane  
Plymouth, MN 55447  
U.S.A.  
Phone: 763.553.1270  
Fax: 763.553.0459  
**1.800.328.3352**



COMFORT



CONFIDENCE



CONTROL