

# **HISTOFLOW**

*Represented By:*

**BCA** *Scientific Inc.*

3395 Brett Road, Mississauga, ON, L5L 4N7

www.BCAscientific.com

Phone:905-820-6125 Fax:905\*820\*8958


## ***VENTED NECROPSY BENCH - 4FT - MODEL # VNB4-1***



This 4ft wide all stainless steel vented necropsy workstation is specially designed for necropsy work in a sitting position.

Positive control of smells and odor is achieved by a downdraft perimeter around the worksurface. The worksurface features a marine edge for containing spills and is removable for cleaning. Adjustable legs allow for the entire workstation to be fixed at a wide range of heights.

Other features include a GFI outlet, fluorescent lighting, lockable instrument drawer, and shelves for your computer, monitor and keyboard.



3-19560 96th Ave. Surrey, B.C.  
Phn:604-888-9056 Fax:604-888-9083

CLIENT	
DESCRIPTION	HISTOFLOW Model VNB4-1 4 Foot Vented Necropsy Bench
DWG NO.	VNB4-1
Sheet	1 of 1

REV	For Review	Date
-----	------------	------