

HISTOFLOW

Represented By:

BCA Scientific Inc.

3395 Brett Road, Mississauga, ON, L5L 4N7

www.BCAscientific.com

Phone:905-820-6125 Fax:905*820*8958


CUSTOM MACROIMAGING TABLE



This all stainless steel macro-imaging workstation is specially designed for imaging work in a sitting position.



The worksurface features a marine edge for containing spills. Adjustable legs allow for the entire workstation to be fixed at a wide range of heights.

Other features include an integral backsplash, lockable instrument drawer and computer shelf.



All Stainless Steel Construction

Scale 1:16



3-19560 96th Ave. Surrey, B.C.
Phn:604-888-9056 Fax:604-888-9083

CLIENT	
DESCRIPTION	Custom Counter Macroimaging Table
DWG NO.	CC38x32AH
Sheet	1 of 1

REV 0

Date