



## GLASS DOOR REFRIGERATORS

Designed to meet the demanding requirements of scientific and laboratory applications.

### Standard Features

- CFC Free Refrigeration with Super Quiet Operation
- Space Saving Design
- Energy Efficient Thermal Glass Doors
- High Density CFC Free Foamed-In-Place Insulation
- Evaporator is Epoxy Coated to Reduce the Potential of Corrosion.
- Large Variety of Sizes (12-72 cu. ft.)
- Bottom Mounted Heavy Duty Compressor
- White Anodized Aluminum Interior
- Stainless Steel Floor
- Self Contained System- No Plumbing Required
- Warranty: 1 Year Parts and Labor;  
5 Year Compressor Parts

### Optional Features

- Casters (Field Installed, Adds 3.5" to Height of Unit)
- Digital Temperature Display
- Electronic Temperature Control with Digital Display
- Audible/ Visual High/Low Alarms and Remote Contacts With Digital Display
- Access Port
- Duplex Outlet
- Chromatography Mast
- Chart Recorder (6", 7 Day)
- Chromatography Models (CH) Which Include:
  - Access Port
  - Duplex Outlet
  - Full Width Shelves



GDM-47




GDM-23

REPRESENTED BY:




3395 Brett Road, Mississauga, ON, L5L 4N7  
www.BCAscientific.com  
Phone: 905-820-6125 / Fax: 905-820-8958



## Hinged Glass Door

Model #	Cu. Ft. Interior	Number of Shelves	Number of Doors	Exterior Dimensions			Crated Weight (lbs.)
				L"	D"	H"	
GDM-12	12	3	1	24 7/8	24 1/2	62 3/8	220
GDM-19	19	4	1	27	25 7/8	75 1/8	265
GDM-23	23	4	1	27	30 7/8	78 5/8	315
GDM-26	26	4	1	30	30 7/8	78 5/8	325
GDM-35	35	8	2	39 1/2	30 7/8	78 5/8	415
GDM-49	49	8	2	54 1/8	30 7/8	78 5/8	470
GDM-72	72	12	3	78 1/8	30 7/8	79 1/4	675

Chromatography Units Include 1 Access Port (specify what side) and 1 Duplex Outlet  
**\*\*Duplex Outlet Has Additional Cord and Plug\*\***



## Sliding Glass Door

Model #	Cu. Ft. Interior	Number of Shelves	Number of Doors	Exterior Dimensions			Crated Weight (lbs.)
				L"	D"	H"	
GDM-33	33	8	2	39 1/2	29 7/8	78 5/8	410
GDM-37	37	8	2	43 1/2	29 7/8	78 5/8	440
GDM-41	41	8	2	47 1/8	29 7/8	78 5/8	440
GDM-45	45	8	2	51 1/8	29 7/8	78 5/8	490
GDM-47	47	8	2	54 1/8	29 7/8	78 5/8	495
GDM-61	61	8	2	70 1/8	29 7/8	78 5/8	585
GDM-69	69	12	3	78 1/8	29 7/8	79 1/4	650

Chromatography Units Include 1 Access Port (specify what side) and 1 Duplex Outlet  
**\*\*Duplex Outlet Has Additional Cord and Plug\*\***



REPRESENTED BY:

**BCASCIENTIFIC**  
 INCORPORATED  
 3395 Brett Road, Mississauga, ON, L5L 4N7  
 www.BCAScientific.com  
 Phone:905-820-6125 / Fax:905-820-8958