

Rees Scientific is proud to introduce our New 5G WiFi system



Our 5G WiFi offers the same great reliability, security and flexibility as the Z3 Wireless System. It's designed to reduce the overall installation process and is great for remote sites. It offers enhanced customer convenience and works with the full line of Rees Scientific sensors.



The 5G spectrum offers improved radio communication and less interference



Has the ability to operate on 2.4 and 5 GHz

Min 2.0°C
Max 8.0°C

Modules with displays show the daily Min and Max in compliance with vaccine storage requirements



Provides a significantly longer and more reliable battery life under most wifi conditions



It offers a higher level of security with a FIPS 140-2 Level 1 validated radio



Offers dependable and consistent data buffering



Has a superior ability to stay online even in the most difficult of conditions



Enterprise WiFi capable option (PEAP-MS-CHAPv2)

Rees Scientific

Continuous Temperature Monitoring



More security and flexibility with the same great service!



MONITORING



PROTECT VALUABLE ASSETS




COMPLIANCE



NOTIFICATION



SECURED WIFI



REPRESENTED BY:



BCASCIENTIFIC INCORPORATED

3395 Brett Road, Mississauga, ON, L5L 4N7

www.BCAScientific.com

Phone:905-820-6125 / Fax:905*820*8958