

Spire Integrated Solutions manufactures the CSI Jewett brand necropsy grossing workstations.

Customizable
Necropsy Grossing Workstation
AR943

This durable and economical workstation, in 3 sizes, with full-width air vent, may be customized to meet the needs of your facility.



MODEL	Dimensions W x D x H inches
AR943	48 x 30 x 41
AR943-60	60 x 30 x 41
AR943-72	72 x 30 x 41

Features

- Made in the USA
- Drain trough includes flushing system and removable dissection grid plates.
- Unit includes a large capacity, multi-directional adjustable ventilation grill.
- Chrome-plated brass faucet has hot and cold water flow control wrist blades & a swivel goose neck spout with integrated vacuum breaker.
- A lever handle provides waste drain control for the sink (deleted with disposer).
- UL listed
- Backdraft ventilation with full-width air vent to connect to the building system
- Large radius inside corners that reduces accumulation of tissue and debris.

Construction

- The sink top, basin, drain board, and backsplash are constructed of welded, 14-gauge stainless steel.
- Heavy-gauge polished stainless steel legs and welded cross-rails constructed to handle heavy use.
- Top, back, sides and shelves are made of 14-gauge, Type 304 stainless steel, polished on exposed and visible surfaces.

Options

- Elevating height adjustment
- Spray hose with palm operated water control
- Task light with switch - various models available
- Base cabinet
- Flushing system
- Factory installed waste disposer with automatic water control
1 hp motors are available (others on request)
- Air, gas and vacuum valves (SR2147-02)
- GFCI outlet
- Scale bracket
- Emergency eye wash
- Customizable



Under work surface
rinse bar.



(Optional) LED magnifying
task light

Drawings

