

Spire Integrated Solutions manufactures the CSI Jewett brand necropsy grossing workstations.

## Necropsy Grossing Workstation JD Series

This compact, shelved workstation, with an optional elevating work surface & adjustable height, is customizable to meet the needs of your facility.



Downdraft ventilation work surface

MODEL	Dimensions W x D x H   inches
JD	48 x 38 x 84

### Features

- Made in the USA
- Removable stainless steel dissecting grid plates located over drain trough.
- Unit includes a large capacity, multi-directional adjustable ventilation grill.
- Exhaust connection with adjustable damper
- Polished surface for easy cleaning
- Wrist-operated, chrome-plated, brass faucets with hot and cold water control valves and swivel goose neck spout.
- Lever-handle, waste drain control for the sink



Perforated work surface with down draft ventilation.

### Construction

- Welded sanitary type 316 stainless steel construction with large radius corners for easy cleaning.
- Heavy-gauge polished stainless steel legs and welded cross-rails constructed to handle heavy use.
- Sink top, basin, drain board, & back splash are made of one-piece, welded, 14-gauge stainless steel.

### Options

- Elevating work surface
- Spray hose with palm operated water control
- Task light with switch - various models available
- Waste disposers - 1/2 to 1.5 hp
- Base cabinet
- Flushing system
- Adjustable height
- Air, gas and vacuum valves (SR2147-02)
- Emergency eye wash
- Customization
- GFCI outlet
- Scale bracket



Large radius corners on interior surfaces for easy cleaning & reduction of tissue & debris accumulation.

Drawings

