

Spire Integrated Solutions manufactures the CSI Jewett brand necropsy workstations.



Fixed workstation with optional elevating workstation

Necropsy Grossing Workstation JNT Series

This popular, durable & economical, fixed workstation with a downdraft work surface is customizable and preferred by histologists nationwide.

DIMENSIONS	Widths	Lengths	Height
JNT Series	30", 36", 42"	48", 60", 72"	Fixed 36" work surface

Features

- Removable stainless steel dissecting grid plates, located over the drain trough, with flushing system.
- Wrist-operated, chrome-plated brass faucet with hot and cold water control valves and swivel goose neck spout.
- Exhaust connection with adjustable damper
- Polished surface for easy cleaning
- Lever-handle, sink waste drain control

Construction

- Sink top, basin, drain board, & back splash are made of one-piece, welded, 14-gauge, type 316 stainless steel.
- Heavy-gauge polished stainless steel legs and welded cross rails to handle heavy usage.
- Welded sanitary type 304 stainless steel construction with large radius inside corners for easy cleaning.

Options

- 316 Stainless steel construction
- Elevating work surface
- Spray hose with palm operated water control
- Waste disposers—1/2 to 1.5 hp
- Task light with switch
- Base cabinet
- Waste disposer
- GFCI Outlet
- Scale bracket
- Emergency eye wash
- Air, gas and vacuum valves (SR2147-02)
- Customizable



Sink with radiused corners with wrist-operated faucet with hot & cold water controls.



Front panel features a coiled spray hose, hot & cold water control knobs, with optional GFCI outlet, & a disposer switch.

Exhaust connection with adjustable damper.



Drawings

