

Cage & Bottle Washers

Model 320LX | 340LX

Hydro-spray washers designed for ease of use, minimal maintenance, and highest cleanability of caging, bottles, debris pans, and other items.



Model	Standard Chamber Size	Standard Unit Size
320LX	48" W x 32" H x 34" D	72" W x 82" H x 40" D
340LX	48" W x 32" H x 48" D	72" W x 82" H x 58" D

Features

- Drop-down Door
- Oscillating Spray Delivery System
- Personnel Safety Features
- PLC Control System with Full Color HMI Touch Screen
- 12 Configurable Cycles
- Neutralization System
- Automatic Self-cleaning Screen
- Temperature Guarantee
- Insulated Construction
- Incoming Utility Gauges
- Interior Illumination
- SS Recirculation Piping and Components
- Non-proprietary parts

Construction

The LYNX brand 320LX and 340LX Cage and Bottle Washers are provided with an all stainless steel drop-down door with a double pane glass window.

Recirculated water is heated by a stainless steel steam coil to achieve guaranteed treatment temperatures. Unit is provided with a stainless steel 7½ HP treatment pump with a direct read pump pressure gauge.

Options

- Pass-thru Unit
- Right or Left Hand Services
- LYNX OPTI-Wash™ Water Saving System
- High Capacity Center Header System
- House Hot Water Heat Exchanger
- Electric Heat System
- LYNX Smart-Cool™ System
- Drain Discharge Cooling - Injected
- Effluent Surge Tank
- Automatic Damper
- Exhaust Fan
- Low Steam Pressure Package
- Non-vented Vapor Condenser
- Pass-thru Door Interlock System
- Air Compressor
- Integral Data Collection System
- Printer
- Seismic Tie Down
- Detergent Dispensing Systems
- pH Neutralization Monitoring System
- Direct DI or pure water final rinse
- Barrier Wall Flange(s)

Accessories

- Cage Wash Rack
- Bottle Wash Rack
- Transfer Cart

Utilities

Electrical (E)	3-60-208/480V, 9HP
Hot Water (HW)	1" NPT, 35 PSI, 30 GPM @ 140° F ± 5° F
Cold Water (CW)*	1" NPT, 50 PSI, 40 GPM @ 50° F ± 5° F
Steam (S)	1½" NPT, 45-80 PSI, 350-600 #/Hr Flow Rate 1" NPT ½"
Condensate (C)	NPT, 80-100 PSI
Air (A)	2-4 CFM Flow Rate
Drain (D)	4" Floor Drain (<i>minimum</i>), 90 GPM Flow Rate 8"
Vent (V)	ID Duct, 200 CFM, Saturated Vapor @ 190° F

*CW - Required if drain discharge cool-down package ordered

Drawings

