

Product Name:	PF021WWW/0M
Description:	<p>White Diamond Series Countertop Freezer Freestanding 1.7 Cu. Ft. Capacity One solid swing door, right hinged* Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in glycol bottle Adjustable operating temperature range: -15°C – -25°C 2 adjustable shelves with guard rail to prevent product from being pushed against back wall Bottom footed shelf to allow air circulation Keyed door lock Magnetic door gasket for positive seal Leveling legs LED interior light Integrated door handle Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Three year parts and labor warranty, plus an additional two year compressor parts only warranty Pyxis[®], Omnicell[®] and AcuDose RX[®] compatible Exterior dimensions: 23 3/4" W x 24" D x 20 3/8" H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Manual defrost 115V, 60 Hz, 1 Amp, 1/8 HP ENERGY STAR CERTIFIED UL/C-UL listed Shipping weight: 106 lbs. *This product is available left hinged, product number PF021WWWLH/0M</p>