

Product Name:	PR051WWWADALH/0
Description:	<p>White Diamond Series Undercounter Refrigerator Built-In-LEFT  HINGED  4.6 Cu. Ft. Capacity  One solid swing door, left hinged  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  2 probes; 1 in air and 1 in glycol bottle  Adjustable operating temperature range: 1°C – 10°C  Three shelves (two adjustable/one fixed) with guard rail on back to prevent product from pushing against back wall  Back wall evaporator cover  Forced draft circulation  Keyed door lock  Magnetic door gasket for positive seal  Heavy duty door handle  Leveling legs  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Three year parts and labor warranty  Pyxis<sup>®</sup>, Omnicell<sup>®</sup> and AcuDose RX<sup>®</sup> compatible  Exterior dimensions: 23 ¾"W x 25 ½"D x 31 15/16"H  Cycle defrost  115V, 60 Hz, 1.74 Amps, 1/8 HP  ENERGY STAR CERTIFIED  ADA compliant and UL/C-UL Listed  Shipping weight: 150 lbs.</p>