

Product Name:	PRF092WWWLH/0
Description:	<p>White Diamond Series Refrigerator & Freezer Combination-LEFT HINGED</p> <p>9 Cu. Ft. capacity (4.6 Ref./4.2 Freezer)</p> <p>Two exterior solid swing doors, left hinged</p> <p>Microprocessor temperature controller, refrigerator and freezer</p> <p>Digital temperature display with sample simulated probe, refrigerator and freezer</p> <p>Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C – -25°C</p> <p>Audible and visual high and low temperature alarms, refrigerator and freezer</p> <p>Remote alarm contacts, refrigerator and freezer</p> <p>2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall</p> <p>Back wall evaporator cover, Refrigerator</p> <p>Forced draft circulation Refrigerator</p> <p>Keyed door locks</p> <p>Magnetic door gasket for positive seal</p> <p>Recessed door handle, freezer; heavy duty door handle, refrigerator</p> <p>Leveling legs</p> <p>Stacking kit</p> <p>Two probe access ports (3/8"); one, refrigerator/one, freezer</p> <p>White powder coated exterior</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Exterior dimensions: 23 3/4"W x 24 1/2"D x 66 3/4"H</p> <p>Manual defrost freezer / cycle defrost refrigerator</p> <p>115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer</p> <p>UL/C-UL listed</p> <p>Shipping weight: 269 lbs.</p> <p>*Please note this unit requires 2 outlets.</p>