

## Washer Disinfecter | Cart & Utensil Washer

**4 sizes to choose from:** SCIENTEK brand Cart and Utensil Washers are heavy-duty, high-capacity hydro-spray washers for sterile processing, emphasizing ease of use, minimal maintenance, and "green" flexibility. Featuring a simple design and use non-proprietary parts.



Model	Standard Chamber Size	Standard Unit Size	Pump
SW6100	41" W x 84" H x 80" D	55.5" W x 91" H x 84" D	7.5 HP
SW6400	48" W x 88" H x 90" D	62.75" W x 95" H x 96" D	10 HP
SW6400L	48" W x 88" H x 128" D	62.75" W x 95" H x 132" D	15 HP
SW6400WL	68.5" W x 88" H x 128" D	83.25" W x 95" H x 132" D	20 HP

### Features

- Heavy-duty Hinge Door
- Modular Construction
- Knock-down Shipment
- Rotating Spray Delivery System
- Personnel Safety Features
- PLC Control System with Full Color HMI Touch Screen
- 10 Configurable Cycles
- Automatic Self-cleaning Screen
- Temperature Guarantee
- Insulated Construction
- Incoming Utility Gauges
- Interior Illumination
- SS Recirculation Piping and Components
- Non-proprietary parts

### Construction

- Heavy-duty stainless steel doors with port view glass windows.
- Recirculated water is heated by a stainless steel steam coil guaranteeing treatment temperatures.
- A color touch screen control panel provides a user-friendly interface for selecting, programming, and monitoring cycles.
- Units are provided with stainless steel treatment pumps.

### Options

- Shallow Base Unit (6")
- Pit or Floor Mounted
- Pass-thru Unit
- Right or Left Hand Services
- House Hot Water Heat Exchanger
- Integral Data Collection System
- Drain Discharge Cooling
- Automatic Damper
- SS Exhaust Fan
- Low Steam Pressure Package
- Pass-thru Door Interlock System
- Air Compressor
- Printer
- Seismic Tie Down
- Detergent Dispensing Systems
- pH Neutralization System
- Split Base for Special Entry
- Utility Enclosure Panels
- Inclined Floor Grating
- Barrier Wall Flange(s)
- Devapourmatic System
- Chemical Level Monitoring System
- Direct quick connect cart washing spray option
- Electric Heating In Place of steam

### Accessories

- Internal Drop in Tilt Ramps
- Exterior Load Ramps (*for floor mounted models*)
- Multi-tier Utensil Wash Cart
- Flat Deck Washing System
- Bottle Washing Spindle Cart

