



ELEVATED COLD STORAGE SOLUTIONS

Space saving design for compact storage needs



Product Portfolio

Extensive product offering ranging from laboratory refrigerators and freezers, to specialized products such as controlled room temperature cabinets for ambient medication storage, and safety-focused solutions for flammable material storage and hazardous locations.



Space Saving Design

Maximize facility space with compact refrigeration, optimizing efficiency and streamlining workflow. Choose from various designs, such as stacked refrigerator and freezer combinations, built-in undercounters, and countertop units - all delivering precise cooling while conserving valuable workspace.



Eco-Friendly Cooling

All refrigerators and freezers utilize natural, hydrocarbon refrigerants that significantly reduce global-warming potential and comply with EPA/SNAP mandates. Choose from over 20 models that meet or exceed Energy Star Certification standards for eco-friendly cooling solutions that prioritize environmental sustainability without compromising performance.

Why Choose Us?

With a focus on precision, reliability, and energy efficiency, Corepoint® Scientific's space-saving refrigerators and freezers are designed to meet the demanding requirements for both scientific research and clinical applications. Maximize efficiency and minimize your footprint with 60+ energy efficient refrigerator and freezer options.

Refrigerators and freezers provide energy efficient, convenient, safe, and reliable performance for optimal storage temperature environments necessary for a wide range of life science, pharmacy, biological, medical, clinical, vaccine and industrial applications.

- Digital temperature displays
- Integrated alarm systems
- Probe access
- Freestanding or Built-in Configurations Available
- Manual, Auto or Controlled Auto Defrost Freezer options



Part Number	Ft' (Ref/Frz)	Application	Setpoint (Ref/Frz)	Defrost (Ref/Frz)	Door (Ref/Frz)	# of Shelves (Ref/Frz)	Exterior Dimensions (w" x d" x h")
Freestanding Refrigerators							
*PR011WWW/O	1	Laboratory	4°C	Cycle	1 Solid	2 Adjustable	17 ¼ x 19 ¼ x 21 ¼
*PR011WWG/O	1	Laboratory	4°C	Cycle	1 Glass	2 Adjustable	17 ¼ x 19 ¼ x 21 ¼
PR031WWW/O	2.5	Laboratory	4°C	Cycle	1 Solid	4 Adjustable	17 ¾ x 19 ½ x 28 ⅞
PR031WWG/O	2.5	Laboratory	4°C	Cycle	1 Glass	4 Adjustable	17 ¾ x 19 ½ x 28 ⅞
FR051WWW/O	5	Flammable Material	4°C	Cycle	1 Solid	3 Adjustable	23 ⅞ x 24 ¼ x 32 ⅜
ER051WWW/O	5	Hazardous Location	4°C	Cycle	1 Solid	3 Adjustable	23 ⅞ x 24 ¼ x 32 ⅜
PR061WWW/O	5.2	Laboratory	4°C	Cycle	1 Solid	3 Adjustable	23 ¾ x 24 x 32 ½
PR061WWG/O	5.2	Laboratory	4°C	Cycle	1 Glass	3 Adjustable	23 ¾ x 24 x 32 ½
PR101WWW/O	10.5	Laboratory	4°C	Cycle	1 Solid	6 Adjustable	20 ¾ x 26 ½ x 59
PR101WWG/O	10.5	Laboratory	4°C	Cycle	1 Glass	6 Adjustable	20 ¾ x 26 ½ x 59
Freestanding Controlled Room Temperature Cabinets							
*CRTPR011WWW/O	1	Pharmacy/Medical	22°C	Cycle	1 Solid	2 Adjustable	17 ¼ x 19 ¼ x 21 ¼
*CRTPR011WWG/O	1	Pharmacy/Medical	22°C	Cycle	1 Glass	2 Adjustable	17 ¼ x 19 ¼ x 21 ¼
CRTPR031WWW/O	2.5	Pharmacy/Medical	22°C	Cycle	1 Solid	4 Adjustable	17 ¾ x 19 ½ x 28 ⅞
CRTPR031WWG/O	2.5	Pharmacy/Medical	22°C	Cycle	1 Glass	4 Adjustable	17 ¾ x 19 ½ x 28 ⅞
CRTPR061WWW/O	5.2	Pharmacy/Medical	22°C	Cycle	1 Solid	3 Adjustable	23 ¾ x 24 x 32 ½
CRTPR061WWG/O	5.2	Pharmacy/Medical	22°C	Cycle	1 Glass	3 Adjustable	23 ¾ x 24 x 32 ½
CRTPR101WWW/O	10.5	Pharmacy/Medical	22°C	Cycle	1 Solid	6 Adjustable	20 ¾ x 26 ½ x 59
CRTPR101WWG/O	10.5	Pharmacy/Medical	22°C	Cycle	1 Glass	6 Adjustable	20 ¾ x 26 ½ x 59
Freestanding Freezers							
PF021WWW/OCAD	1	Pharmacy/Medical	-20°C	Controlled Auto	1 Solid	2 Adjustable	23 ¾ x 24 x 20 ⅜
*PF021WWW/OA	1.3	Laboratory	-20°C	Auto	1 Solid	2 Adjustable	23 ¾ x 24 x 20 ⅜
*PF021WWW/OM	1.7	Laboratory	-20°C	Manual	1 Solid	2 Adjustable	23 ¾ x 24 x 20 ⅜
PF051WWW/OM20	4	Laboratory	-20°C	Manual	1 Solid	2 Fixed+3 Baskets	23 ⅞ x 24 ¼ x 32 ⅜
PF051WWW/OM30	4	Laboratory	-30°C	Manual	1 Solid	3 Adjustable	23 ¾ x 24 x 32 ½
PF051WWW/OM40	4	Laboratory	-40°C	Manual	1 Solid	3 Adjustable	23 ¾ x 24 x 32 ½
FF051WWW/OM	4	Flammable Material	-20°C	Manual	1 Solid	2 Fixed+3 Baskets	23 ⅞ x 24 ¼ x 32 ⅜
EF051WWW/OM	4	Hazardous Location	-20°C	Manual	1 Solid	2 Fixed+3 Baskets	23 ⅞ x 24 ¼ x 32 ⅜
Freestanding Combinations							
PRF072WWG/OCAD	5.2/1	Pharmacy/Medical	4°C/-20°C	Cycle/Controlled Auto	1 Glass/1 Solid	3 Adjustable/2 Adjustable	23 ¾ x 24 ¼ x 53 ⅞
PRF072WWW/OCAD	5.2/1	Pharmacy/Medical	4°C/-20°C	Cycle/Controlled Auto	2 Solid	3 Adjustable/2 Adjustable	23 ¾ x 24 x 53 ⅞
PRF072WWG/OA	5.2/1.3	Laboratory	4°C/-20°C	Cycle/Auto	1 Glass/1 Solid	3 Adjustable/2 Adjustable	23 ¾ x 24 x 53 ⅞
PRF072WWW/OA	5.2/1.3	Laboratory	4°C/-20°C	Cycle/Auto	2 Solid	3 Adjustable/2 Adjustable	23 ¾ x 24 x 53 ⅞
PRF072WWG/OM	5.2/1.7	Laboratory	4°C/-20°C	Cycle/Manual	1 Glass/1 Solid	3 Adjustable/2 Adjustable	23 ¾ x 24 x 53 ⅞
PRF072WWW/OM	5.2/1.7	Laboratory	4°C/-20°C	Cycle/Manual	2 Solid	3 Adjustable/2 Adjustable	23 ¾ x 24 x 53 ⅞
PRFI02WWG/O	5.2/4	Laboratory	4°C/-20°C	Cycle/Manual	1 Glass/1 Solid	3 Adjustable/2 Fixed+3 Baskets	23 ¾ x 24 x 65
PRFI02WWW/O	5.2/4	Laboratory	4°C/-20°C	Cycle/Manual	2 Solid	3 Adjustable/2 Fixed+3 Baskets	23 ¾ x 24 x 65
PRFI02WWG/O30	5.2/4	Laboratory	4°C/-30°C	Cycle/Manual	1 Glass/1 Solid	3 Adjustable/3 Adjustable	23 ¾ x 24 x 65
PRFI02WWW/O30	5.2/4	Laboratory	4°C/-30°C	Cycle/Manual	2 Solid	3 Adjustable/3 Adjustable	23 ¾ x 24 x 65
PRFI02WWG/O40	5.2/4	Laboratory	4°C/-40°C	Cycle/Manual	1 Glass/1 Solid	3 Adjustable/3 Adjustable	23 ¾ x 24 x 65
PRFI02WWW/O40	5.2/4	Laboratory	4°C/-40°C	Cycle/Manual	2 Solid	3 Adjustable/3 Adjustable	23 ¾ x 24 x 65
PRFI22WWG/OCAD	10.5/1	Pharmacy/Medical	4°C/-20°C	Cycle/Controlled Auto	1 Glass/1 Solid	6 Adjustable/2 Adjustable	23 ¾ x 26 ½ x 80 ¼
PRFI22WWW/OCAD	10.5/1	Pharmacy/Medical	4°C/-20°C	Cycle/Controlled Auto	2 Solid	6 Adjustable/2 Adjustable	23 ¾ x 26 ½ x 80 ¼
PRFI22WWG/OA	10.5/1.3	Laboratory	4°C/-20°C	Cycle/Auto	1 Glass/1 Solid	6 Adjustable/2 Adjustable	23 ¾ x 26 ½ x 80 ¼
PRFI22WWW/OA	10.5/1.3	Laboratory	4°C/-20°C	Cycle/Auto	2 Solid	6 Adjustable/2 Adjustable	23 ¾ x 26 ½ x 80 ¼
PRFI22WWG/OM	10.5/1.7	Laboratory	4°C/-20°C	Cycle/Manual	1 Glass/1 Solid	6 Adjustable/2 Adjustable	23 ¾ x 26 ½ x 80 ¼
PRFI22WWW/OM	10.5/1.7	Laboratory	4°C/-20°C	Cycle/Manual	2 Solid	6 Adjustable/2 Adjustable	23 ¾ x 26 ½ x 80 ¼
Built-in Refrigerators							
*PR031WWW/OFB	2.5	Laboratory	4°C	Cycle	1 Solid	4 Adjustable	17 ¾ x 22 ¾ x 30 ¾
*PR051SSS/O	4.6	Laboratory	4°C	Cycle	1 Solid	2 Adjustable/1 Fixed	23 ¾ x 24 ½ x 33 ⅞
*PR051WWW/O	4.6	Laboratory	4°C	Cycle	1 Solid	2 Adjustable/1 Fixed	23 ¾ x 24 ½ x 33 ⅞
*PR051WWG/O	4.6	Laboratory	4°C	Cycle	1 Glass	2 Adjustable/1 Fixed	23 ¾ x 25 ½ x 33 ⅞
*PR051WWWADA/O	4.6	Laboratory	4°C	Cycle	1 Solid	2 Adjustable/1 Fixed	23 ¾ x 24 ½ x 31 15/16
*PR051WWGADA/O	4.6	Laboratory	4°C	Cycle	1 Glass	2 Adjustable/1 Fixed	23 ¾ x 25 ½ x 31 15/16
Built-in Controlled Room Temperature Cabinets							
*CRTPR031WWW/OFB	2.5	Pharmacy/Medical	22°C	Cycle	1 Solid	4 Adjustable	17 ¾ x 22 ¾ x 30 ¾
*CRTPR051WWW/O	4.6	Pharmacy/Medical	22°C	Cycle	1 Solid	2 Adjustable/1 Fixed	23 ¾ x 24 ½ x 33 ⅞
*CRTPR051WWG/O	4.6	Pharmacy/Medical	22°C	Cycle	1 Glass	2 Adjustable/1 Fixed	23 ¾ x 25 ½ x 33 ⅞
Built-in Freezers							
*PF051WWW/OM	4.2	Laboratory	-20°C	Manual	1 Solid	2 Fixed	23 ¾ x 24 ½ x 33 ⅞
*PF051WWWADA/OM	4.2	Laboratory	-20°C	Manual	1 Solid	2 Fixed	23 ¾ x 24 ½ x 31 15/16
*PF051SSS/OM	4.2	Laboratory	-20°C	Manual	1 Solid	2 Fixed	23 ¾ x 24 ½ x 33 ⅞
*PF051WWW/OA	4.2	Laboratory	-20°C	Auto	1 Solid	3 Adjustable	23 ¾ x 24 ½ x 33 ⅞
Built-in Combinations							
*PRF092WWG/O	4.6/4.2	Laboratory	4°C/-20°C	Cycle/Manual	1 Glass/1 Solid	2 Adjustable, 1 Fixed/2 Fixed	23 ¾ x 25 ½ x 66 ¾
*PRF092WWW/O	4.6/4.2	Laboratory	4°C/-20°C	Cycle/Manual	2 Solid	2 Adjustable, 1 Fixed/2 Fixed	23 ¾ x 24 ½ x 66 ¾
*PRF092SSS/O	4.5/4.2	Laboratory	4°C/-20°C	Cycle/Manual	2 Solid	2 Adjustable, 1 Fixed/2 Fixed	23 ¾ x 24 ½ x 66 ¾

*Left hinge model available