



Spire Integrated Solutions manufactures the CSI Jewett brand of healthcare equipment.

CWC Series Warming Cabinets

From stationary cabinets to mobile units, the CSI Jewett brand provides premium warming technology designed for both blankets and fluids.



CWC-15 model warming cabinet with cart.

Model	Exterior Cabinet Dimensions W x D x H inches	Interior Cabinet Dimensions W x D x H inches	Weight lbs.	Interior Volume cu. ft	Glass Door Type	Standing or Countertop
CWC-15	30 x 24 x 28	25.5 x 17.5 x 15.5	152	4' ³	Solid or dual pane, low E	Both
CWC-28	20 x 24 x 41	25.5 x 17.5 x 28.5	152	7.4' ³	Solid or dual pane, low E	Both
CWC-41	30 x 24 x 54	25.5 x 17.5 x 41.5	263	10.8' ³	Solid or dual pane, low E	Both
CWC-61	30 x 24 x 74	25.5 x 17.5 x 61.5	510	15.9' ³	Solid or dual pane, low E	Both w/ opt. casters
CWC-1515	30 x 24 x 48	25.5 x 17.5 x 15.5	310	4.3' ³	Solid or dual pane, low E	Both w/ opt. casters
CWC-1528	30 x 24 x 61	25.5 x 17.5 x 28.5	363	7.4' ³	Solid or dual pane, low E	Both w/ opt. casters
CWC-1541	30 x 24 x 74	25.5 x 17.5 x 41.5	414	10.75' ³	Solid or dual pane, low E	Both w/ opt. casters
CWC-2828	30 x 24 x 74	25.5 x 17.5 x 28.5	422	7.4' ³	Solid or dual pane, low E	Both w/ opt. casters

Features

- Fully programmable user controls for your specific requirements.
- Large digital screen with individual cabinet controls.
- Fully insulated and enclosed walls and doors.
- Adjustable perforated stainless steel shelves prevent bottle tipping.
- Positive door latches with self-closing hinges.
- Continuous forced fan hot air convection over electric heating element.
- Hospital grade power cord.
- Right or left hinge from the factory, field reversible.
- Automatic high temperature cutout and reset.
- Data logger with removable USB memory.
- Leveling feet
- TUV Listed To UL 61010A-1/ & CSA C22.2, No. 61010-1

Construction

- Stainless steel construction for durability and longevity.

Options

- 18" depth casters
- Insulated, tempered safety glass doors
- Electronic door locks
- Keyed or control panel lock out
- Seismic anchors
- Stainless steel trim
- Additional shelves
- Door ajar alarm
- Pass-through configurations
- Cart
- 220 volt electrical



Fully programmable user controls for your specific requirements. (Optional) Control panel lock out



(Optional) Keyed lock out



Adjustable perforated stainless steel shelves prevent bottle tipping.



Warming Cabinet Cart

The CSI Jewett brand of warming cabinets also offers an optional mobile transport and storage cart for use with our CWC series. These Carts feature full stainless steel construction, casters with both swivel and locking features, and offer storage, mobility and flexibility for our smaller model fluid and blanket warmers.

Features

- (2) shelves with access holes to securely mount the cabinet to the cart on the top shelf.
- Bottom shelf available for storing additional blankets or supplies.
- Wheels